



Recombinant *Drosophila melanogaster* RNA-binding protein Musashi homolog Rbp6 (Rbp6)

Product Code	CSB-EP893335DLU-B
Abbreviation	Rbp6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VVE5
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDIKIEQQQQ QQQVELGPCS PSEVPNDPGK MFIGGLSWQT SPESLRDYFG RYGDISEAMV MKDPTTRRSR GFGFVTFSDP NSVDKVLQGG THELDGKKVD PKVAFPRRAH PKMVTRTKKI FVGGLSAPTT LEDVKSYPEQ FGPIEDAMLM FDKQTNRHRG FGFVTFQSED VVDKVCEIHF HEINNKMVEC KKAQPKEVML PANLAKTRAA GRSAYDNIMW GLGTLPEGFP AAAYAAYAAG RGYSGYPSFG LPYPTVDLTN GQLTPNNNTP LLTQALSVMN NYQAAAAAAH SYGPSASPHA AAANTATRQG FPSTNSPGPT IDMYSSAAD NVGYVQATSP QPSGFPIAVS RAPLNYSPIY SLNFAMANI
Source	E.coli
Target Names	Rbp6
Protein Names	Recommended name: RNA-binding protein Musashi homolog Rbp6
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.