



# Recombinant *Drosophila melanogaster* Protein painting of fourth (Pof)

<b>Product Code</b>	CSB-EP893358DLU-B
<b>Abbreviation</b>	Pof
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9W123
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDSKRAALES GDGPDAKRLD TTDDQDKEAS GGDGSQVMLA KHVAPYTGHG CTPPMESYLF EPTPAGSPLL PWKTSVDLDN DAELEKPTSD KKPD TAKLSR RELAKMRREH TLRALALERE LTNKPGQTPA SEVLLVRFDP PEITAPMLAG LSKDIRDVVL PISVAPRYCL VHLKAGADVE ATICDINRVR FGTGHLRAEL KPFSDEEQAE FIDPCSLYVG NIPFNMTTSA IKAYFANAMR VDIGVLKREK RARYAFVRYA SPDQTMFAFK ELVDSPLNSR TLTVRYRRLR KRAGMPMVQC ATSFQALQSP NGDDDNTDCK VISPPPLESI IISDSDNCSD SSGNGKEDGK RKKKINEQER EIEKLKRQMA EYGAIKSLQ FRQNSLEDTF IPDLTPKVEP SVNPTGCLLG SNAVHLMRDI KKECDYLGIP DPVPATKPTT QAQDDSQKKA KRSCFGRLFT GPFRRGTSAM KTADEYEKDD RLEELYAQLR RDPDP
<b>Source</b>	E.coli
<b>Target Names</b>	Pof
<b>Protein Names</b>	Recommended name: Protein painting of fourth Alternative name(s): Zeste-interacting protein 16 Short name= Zip16
<b>Expression Region</b>	1-495
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.