



Recombinant *Drosophila melanogaster* OCIA domain-containing protein 1 (asrij)

Product Code	CSB-EP893362DLU-B
Abbreviation	asrij
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W1X9
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MDSPLNDGSH HPPPHAPHPL ADYQFSAEEV KALRECNTEs FFQRSLPFGT GLGLLAYFGV KNGYLQGHVK YGAVPKVVMG VILGYFVGKF SYQQKCAEKI MRLPNSHLGE LLRQRRQGGG VISSITPDEN LGRAFTLAPF SPSSADVYS EAYQPGRSTS LNLDTESRPT LSGLDDIYRP TLDSGSMLEA ELPLEPSKPG QSYEDLRRRN REEYSKHQQS PYSRPEYPPV AVQQRPEQA QSEPAGRKNQ YGDSWTD
Source	<i>E. coli</i>
Target Names	asrij
Protein Names	Recommended name: OCIA domain-containing protein 1
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.