



Recombinant Drosophila melanogaster Probable U3 small nucleolar RNA-associated protein 11 (CG1789)

Product Code	CSB-BP893365DLU
Abbreviation	CG1789
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W3C0
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MSSWKNASKS NQKVHRERHQ PESRQHLGFL EKKKDYKKRA IDAQKKQKTL KLLYRRAQNK NPDEFYHHMI NSKLSNDEHH EKDTKDEHTP EQLALMQTQD LKYYVMKRTM ERKKIERLKA SLVDVDAIKG AANKRIKFDE DGYMELDMGE WLLKHDDEEE EEEKKKSKT TEPNPIKQQR INELKKREQR ERELGIVQQK IQLKNTLKQP RLLKPRKIKP GTKDSAPVYK FRYERKK
Source	Baculovirus
Target Names	CG1789
Protein Names	Recommended name: Probable U3 small nucleolar RNA-associated protein 11 Short name= U3 snoRNA-associated protein 11
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.