



Recombinant Drosophila melanogaster Transcription initiation factor IIA subunit 2-2 (TfIIA-S-2)

Product Code	CSB-MP893373DLU
Abbreviation	TfIIA-S-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W5B9
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MNYQHRYRATT LGRTLQDTLD EMMERGDITK KIANLVLLRY DKSISTALKD HGTSNMSFTA ERLETFRCCD NVWTLILKDA EFREDQHSLK VDVVKIVACL GTDNGNE
Source	Mammalian cell
Target Names	TfIIA-S-2
Protein Names	Recommended name: Transcription initiation factor IIA subunit 2-2 Alternative name(s): General transcription factor IIA subunit 2-2 Transcription initiation factor IIA gamma-2 chain Short name= TFIIA-gamma-2
Expression Region	1-107
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.