



Recombinant Drosophila melanogaster Histone H3-like centromeric protein cid (cid)

Product Code	CSB-MP893714DLU
Abbreviation	cid
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V6Q2
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MPRHSRAKRA PRPSANNSKS PNDDDTAFRS PEPEDGTDYG LEFTTSQLTL QDNNRRSSTL RRDAGRQPA ARDSSTSGEE EDQENRYPTT RSPQTRRMTV QQESKTRAAG PVAAQNQTRR RKAANPMSRA KRMDREIRRL QHHPGTLIPK LPFSRLVREF IVKYSDDDEPL RVTEGALLAM QESCEMYLTQ RLADSYMLTK HRNRVTLEVR DMALMAYICD RGRQF
Source	Mammalian cell
Target Names	cid
Protein Names	Recommended name: Histone H3-like centromeric protein cid Alternative name(s): CENP-A homolog Centromere identifier protein
Expression Region	1-225
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.