



Recombinant *Drosophila melanogaster* Proteasome subunit alpha type-4-like (Prosalpha3T)

Product Code	CSB-MP893722DLU
Abbreviation	Prosalpha3T
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VA12
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MARFFDSRTT IFSPEGRLYQ VEYAMEAASQ SGTCVGLLAK NGVLLATERS VDKLMDSIP VPRISWLNEN IACCATGNTA DGNVLVNLQR MIAQQYQFNF GEMIPCEQLV TNLCDIKQAY TQYGGKRPFV VSFLYMGWDC RFGFQLYQSD PSGNYSGWKA TCIGRKSGAA MEMLQKELFS KGYVSPSVEE AKDVAIKVMG MTLGRDSLTP EKLEIAFVQR YGNTTVFHIL EKNEIHRLE RNNNLKRRVG S
Source	Mammalian cell
Target Names	Prosalpha3T
Protein Names	Recommended name: Proteasome subunit alpha type-4-like EC= 3.4.25.1
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.