



Recombinant Drosophila melanogaster Probable tRNA (guanine (26)-N (2))-dimethyltransferase (CG6388), partial

| | |
|--------------------------|---|
| Product Code | CSB-BP893768DLU |
| Abbreviation | CG6388 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9VK89 |
| Product Type | Recombinant Protein |
| Immunogen Species | Drosophila melanogaster (Fruit fly) |
| Purity | >85% (SDS-PAGE) |
| Source | Baculovirus |
| Target Names | CG6388 |
| Protein Names | Recommended name: Probable tRNA (guanine(26)-N(2))-dimethyltransferase EC= 2.1.1.216 Alternative name(s): tRNA 2,2-dimethylguanosine-26 methyltransferase tRNA(guanine-26,N(2)-N(2)) methyltransferase tRNA(m(2,2)G26)dimethyltransferas |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |