



Recombinant Drosophila melanogaster Protein AATF-like (CG11188)

Product Code	CSB-EP893776DLU-B
Abbreviation	CG11188
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VM95
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MLRKSkkQPQ TVAEKVSKLL AHPNESDSA E DSDFDVATGP RLVD FEEEE EY DLPDARSTDF RKRNVKLLSE QSDRYK GKIS SRKELDDDED KDDEQEVS YE ESDEDDENLT NFKQKLNAGG AEDSEEETAA GHSESGEESE EIESNL TDFK KKFEAGDFKY DDDEEEDDDS EEEDHSQESE GDDDDSEDDE AEDDAIKPSD VMSKTNHQAE IQKGLAVQNNQ LRIWERLLEL RINTQKFTSK ANQLPAPETL TKLASESDEL QSVLNEAQER SSKLLQQLLS LQSALHQQYS EMKKS VKRKQ PTEDSGPAVK KFGSVLQSNF QQMIGYRNEV LLKWDDRTKL LTPGAGVKRK SLQEDYDIK KIGSALANRE ALVEKSQTPK NSQAEQQENT PVQRLKHIYD DSDFYHQQLR ELIEYKASTS SNMSEITKQF VELQKLRQKM KKKVDTRASK GRKLRYVVHN KLINFMAPNE SSDWTDASKS ELYKSLFV
Source	E.coli
Target Names	Aatf
Protein Names	Recommended name: Protein AATF-like
Expression Region	1-488
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.