



Recombinant *Drosophila melanogaster* Intraflagellar transport protein 57 homolog (CG8853)

Product Code	CSB-BP893790DLU
Abbreviation	CG8853
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VQS5
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MQQDDEQEKS QQLQNFQSDD LLEKLLKLLNY EKHLLEKFKL KPLSRFYFVK ATNPGEQFYM FTLICWWLCK KLGKDMERPQ EYDDPNTVAA NIIKLLGEID VPVDFQPNKL IRGAGPICLS VLEVLSTQAC KVAQVGYQKL HIAQEEEEFLG DYLEDNAEII LEKLEDEQNA AALSDDSDME LEAHNFRQLN WLNRPQKKS GDVNLDERNP ELDARMSDHQ EWHLELERVL PQLKVFVKAD ARDWRTHISQ METLKTNILE NSDTAEAQLK KXHSEFTFDL EKIESREKHL NNELKPLIQQ FKELSIELST IQYAQNQLQE DMEKQTAELN EVMMEQELKK EEMERRGQAM SDGSSVQQIR KAIAKLKDDT AQLNLEVALL VHAHDQDIVV RQLQQNQTTD LANNP
Source	Baculovirus
Target Names	IFT57
Protein Names	Recommended name: Intraflagellar transport protein 57 homolog
Expression Region	1-405
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.