



Recombinant Rat Podocalyxin (Podxl), partial

Product Code	CSB-YP893876RA
Relevance	Involved in the regulation of both adhesion and cell morphology and cancer progression. Function as an anti-adhesive molecule that maintains an open filtration pathway between neighboring foot processes in the podocyte by charge repulsion. Acts as a pro-adhesive molecule, enhancing the adherence of cells to immobilized ligands, increasing the rate of migration and cell-cell contacts in an integrin-dependent manner. Induces the formation of apical actin-dependent microvilli. Involved in the formation of a preapical plasma membrane subdomain to set up initial epithelial polarization and the apical lumen formation during renal tubulogenesis. Plays a role in cancer development and aggressiveness by inducing cell migration and invasion through its interaction with the actin-binding protein EZR. Affects EZR-dependent signaling events, leading to increased activities of the MAPK and PI3K pathways in cancer cells.
Abbreviation	Recombinant Rat Podxl protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WTQ2
Alias	Podocalyxin-like protein 1 ;PC ;PCLP-1
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	QDNGNKTDTSIDTSIDQNQDKPATNQPSNATPKSSVQPPTPTSISTSSPDKA TQSSNSSVTTTTSDSTTDRTSSSTSTVPTTNSGQTVSSGGKSSDKITTALPTTL GPVNASSQPTDLNTSTKLPSTPTTNSTASPHQPVSHSEGQHTTVQSSSASVS SSDNTLLWILTTSKPTGTSEGTQPIAISTPGITTPVSTPLQPTGSPGGTESVPT TEEFTHSTSSWTPVVSQGPSTPSSTWTSGSYKCLKDPAIKPHEELLILNLTRDS FCKGSPNERFLELLCHSAKASFKAEDSCALELAPILDNQAVAVKRIVETKLS PKAVFELLKDKWDDLTEAGVIDIHLGKEGPPEVNEDRFS
Research Area	Others
Source	Yeast
Target Names	Podxl
Protein Names	Recommended name: Podocalyxin Alternative name(s): Podocalyxin-like protein 1 Short name= PC Short name= PCLP-1
Expression Region	25-386aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.




Tag Info N-terminal 6xHis-tagged

Mol. Weight 39.8kDa

Protein Length Extracellular Domain

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.