



Recombinant Mouse Cullin-1 (Cul1)

Product Code	CSB-BP893880MO
Abbreviation	Cul1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WTX6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MSSNRSQNP H GLKQIGLDQI WDDL RAGIQQ VYTRQ SMAKS RYMELYTHVY NYCTSVHQSN QARGAGVPPS KSKKGQTPGG AQFVGLLELYK RLKEFLKNYL TNULLKDGEDL MDESVLKFYT QWEDYRFSS KVLNGICAYL NRHWVRRECD EGRKGIYEIY SLALVTWRDC LFRPLNKQVT NAVLK LIEKE RNGETINTRL ISGVVQSYVE LGLNEDDAFA KGPTLT VYKE SFESQFLADT ERFYTRESTE FLQQNPVTEY MKKAEARLLE EQRRVQVYLH ESTQDELARK CEQV LIEKHL EIFHTEFQNL LDADKNEDLG RMYNLVSR IQ DGLGELK KLL ETHIHNQGLA AIEKCGE AAL NDPKMYVQTV LDVHKKYNAL VMSAFNNDAG FVAALDKACG RFINNNAVTK MAQSSSKSPE LLARYCDSLL KKSSKNPEEA ELEDTLNQVM VVFKYIEDKD VFQKFYAKML AKRLVHQNSA SDDAEASMIS KLKQACGF EY TSKLQRMFQD IGVSKDLNEQ FKKHLTNSEP LDLDFSIQVL SSGSWPFQQS CTFALPSELE RSYQRFTAFY ASRHSGRKL T WLYQLSKGEL VTNC FKNRYT LQASTFQMAI LLQYNTEDAY TVQQLT DSTQ IKMDILAQVL QILLKSKLLV LEDENANVDE VELKPDTLIK LYLGYKNK KL RVNINVPMKT EQKQE QETTH KNIEEDR KLL IQAAIVRIMK MRKVLKHQQL LGEVLTQLSS RFKPRVPVIK KCIDILIEKE YLERVDGEKD TYSYLA
Source	Baculovirus
Target Names	Cul1
Protein Names	Recommended name: Cullin-1 Short name= CUL-1
Expression Region	1-776
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.