



Recombinant Microcin J25-processing protein mcjB (mcjB)

Product Code	CSB-MP893972ENL
Abbreviation	mcjB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9X2V8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	MIRYCLTSYR EDLVILDIIN DSFSIVPDAG SLLKERDKLL KEFPQLSYFF DSEYHIGSVS RNSDTSFLEE RWFLPEPDKT LYKCSLFKRF ILLLKVFYYS WNIEKKGMAW IFISNKKENR LYSLNEEHLI RKEISNLSII FHLNIFKSDC LTYSYALKRI LNSRNIDAHL VIGVRTQPFY SHSWVEVGGQ VINDAPNMRD KLSVIAEI
Source	Mammalian cell
Target Names	mcjB
Protein Names	Recommended name: Microcin J25-processing protein mcjB
Expression Region	1-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.