



Recombinant Arabidopsis thaliana Putative cytokinin riboside 5'-monophosphate phosphoribohydrolase LOG9 (LOG9)

Product Code	CSB-MP894036DOA
Abbreviation	LOG9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XH06
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MHIEHISGET VGEVRIVSDM HERKATMAQE AGAFIALLGE RYETMEELLE MITWAQLGIH KKTVGLLNVD GYYNNLLAFF DTGVEEGFIK QGACNIVVSA PSARELMEM ELYTPSHKYI ASHQSWKVEP LGDYPLLNEN KPQ
Source	Mammalian cell
Target Names	LOG9
Protein Names	Recommended name: Putative cytokinin riboside 5'-monophosphate phosphoribohydrolase LOG9 EC= 3.2.2.n1 Alternative name(s): Protein LONELY GUY 9
Expression Region	1-143
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.