



Recombinant Ribosome biogenesis regulatory protein homolog (rrbs-1)

Product Code	CSB-MP894084CXY
Abbreviation	rrbs-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XVT0
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MADKSTEVEK AIDPIIDLGN LLFIDREPIQ GDASEGLEER ARKNTQLLFN NIWQLEQKRV EEAIIVTLPP ATYRLPREKR LPEKKEPTKW EKYAAEKGIE KRKDKKVF D EATKEWKPTY GYRRGNDDTK DWLIEIPDNA EDPNKDFFAE RREKKKVERVA KNDMQRMKNL ARQMKT TVKS GPSTEKMIGV GVDAKEKSKQ QVRFAVDRAK LATASAGKFQ EGLKGEKANV KTGKKRKFEA NEAPVSGEKE RALQILQRMK SKKAKIVEEK AAVAGPLRE KKEKSERKGA KDQTRQKSQI HRQQWFKNKV DGKKKGTGGA GKKGANKAKA RKG
Source	Mammalian cell
Target Names	rrbs-1
Protein Names	Recommended name: Ribosome biogenesis regulatory protein homolog
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.