



Recombinant Human Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit beta isoform (PIK3CB), partial

Product Code	CSB-EP017999HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P42338
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MCFSFIMPPA MADILDIWAV DSQIASDGS I PVDFLLPTGI YIQLEVPREA TISYIKQMLW KQVHNYPMFN LLM DIDSYMF ACVNQTAVYE ELEDETRRLC DVRPFLPVLK LVTRSCDPGE KLDSKIGVLI GKGLHEFDSL KDPEVNEFRR KMRKFSEEKI LSLVGLSWMD WLKQTYPPEH EPSIPENLED KLYGGKLIVA VHFENCQDVF SFQVSPNMNP IKVNELAIQK RLTIHGKEDE VSPYDYVLQV SGRVEYVFGD HPLIQFYIR NCVMNRALPH FILVECCKIK KMYEQEMIAI EAAINRNSSN LPLPLPPKKT RIISHVWENN NPFQIVLVKG NKLNTEETVK VHVRAGLFHG TELLCKTIVS SEVSGKNDHI WNEPLEFDIN ICDLPRMARL CFAVYAVLDK VKTKKSTKI NPSKYQTIRK AGKVHYPVAW VNTMVFDFKG QLRTGDIILH SWSSFPDELE EMLNPMGT VQ TNPYTENATA LHVKFPENKK QPYYP PFDK IIEKAAEIAS SDSANVSSRG GKKFLPVLKE ILDRDPLSQL CENEMDLIWT LRQDCREIFP QSLPKLLLSI KWNKLEDVAQ LQALLQIWPK LPPREALELL DFNYPDQYVR EYAVGCLRQM SDEELSQYLL QLVQVLKYEP FLDCALSRFL LERALGNRRI GQFLFWHLRS EVHIPAVSVQ FGVILEAYCR GSVGHMKVLS KQVEALNKLK TLNSLIKLNA VKLNRAKGKE AMHTCLKQSA YREALSDLQS PLNPCVILSE LYVEKCKYMD SKMKPLWL VY NNVKVFGE DSV GVIFKNGDDL RQDMLTLQML RLMDLLWKEA GLDLRMLPYG CLATGDRSGL IEVVSTSETI ADIQLNSSNV AAAAAFNKDA LLNWLKEYNS GDDL DRAIEE FTLSCAGYCV ASYVLGIGDR HSDNIMVKKT GQLFHIDFGH ILGNFKSKFG IKRERVPIFIL TYDFIHVIQQ GKTGNTEKFG RFRQCCEDAY LILRRHGNLF ITL FALMLTA GLPELTSVKD IQY LKDSLAL GKSEEEALKQ FKQKFDEALR ESWTTKVNWM AHTVRKDYRS</p>
Source	E.coli
Target Names	PIK3CB
Protein Names	Recommended name: Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit beta isoform Short name= PI3-kinase subunit beta Short name= PI3K-beta Short name= PI3Kbeta Short name= PtdIns-3-kinase subunit beta EC=
Expression Region	1-1070
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.



Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.