



Recombinant Drosophila melanogaster Protein lethal (2)k10201 (I (2)k10201)

Product Code	CSB-MP894182DLU
Abbreviation	I(2)k10201
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V574
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDSEAAGVAL LSREQIIRIL SEVPAGFIKP TDPPAPSLAP FKKLGVLVDI DDILGEQDAA AARIRDNIDK EQSYSCVECR KMLPTAHL LD LHITEQHDCY FAASVERGDK PMFSCFLEEC TIKFHTARQR KDHCIIITHKL PANYRFDHSK NRGKQKHQ GK SKPNSMEVDE VIEETKSLPY VKAFSFGHQ T QRSFYTGKDQ RSGKTLDDVQ AMKEAINDIL D
Source	Mammalian cell
Target Names	I(2)k10201
Protein Names	Recommended name: Protein lethal(2)k10201 Alternative name(s): Wunen region B protein
Expression Region	1-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.