



# Recombinant *Drosophila melanogaster* GDP-fucose protein O-fucosyltransferase 1 (O-fut1)

<b>Product Code</b>	CSB-BP894186DLU
<b>Abbreviation</b>	O-fut1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9V6X7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QLGGDPN GYLTYCPCMG RFGNQADHFL GSLAFKALN RTLILPPWVE YRRGELRSRQ VPFNTYFEVE PLKEYHRVIT MADFMWHLAD DIWPESERVS FCYKERYSLQ QEKNPDKPN CHAKDGNPFG PFWDTFHIDF VRSEFYAPLH FDVHHSNEAA KWQTKYPAES YPVLAFTGAP ASFPVQLENC KLQRYLQWSQ RYREASKDFI REQLPRGAFL GIHLRNGIDW VRACEHVKDS QHLFASPQCL GYKNERGALY PELC MPSKEA IIRQLKRTIK NVRQTQPDNE IKS FV SASDS NHMIGELNTA LSRMGISVHK LPEDDPYLDL AILGQSNHFI GNCISSYSAF EKRRERDVHGF PSYFWGFPKE KDRKHTNVHE EL
<b>Source</b>	Baculovirus
<b>Target Names</b>	O-fut1
<b>Protein Names</b>	Recommended name: GDP-fucose protein O-fucosyltransferase 1 EC= 2.4.1.221 Alternative name(s): Neurotic protein Peptide-O-fucosyltransferase 1 Short name= O-FucT-1
<b>Expression Region</b>	24-402
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.