



Recombinant *Drosophila melanogaster* COP9 signalosome complex subunit 6 (CSN6)

Product Code	CSB-MP894206DLU
Abbreviation	CSN6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VCY3
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MEQMEVDVDM SAKPSTSSSA AAGSSMAVDK TADQNPQPQG NIMAAAGTSG SVTISLHPLV IMNISEHWTR FRAQHGEPRQ VYGALIGKQK GRNIEIMNSF ELKTDVIGDE TVINKDYINK KEQQYKQVFS DLDFIGWYTT GDNPTADDIK IQRQIAAINE CPIMLQLNPL SRSVDHLPLK LFESLIDLVD GEATMLFVPL TYTLATEEAE RIGVDHVARM TSNESGEKSV VAEHLVAQDS AIKMLNTRIK IVLQYIRDVE AGKLRANQEI LREAYALCHR LPVMQVPAFQ EEFYTQCNDV GLISYLGTLT KGCNDMHHFV NKFNMPLYDRQ GSARRMRGLY Y
Source	Mammalian cell
Target Names	CSN6
Protein Names	Recommended name: COP9 signalosome complex subunit 6 Short name= Dch6 Short name= Signalosome subunit 6
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.