



Recombinant *Drosophila melanogaster* Tissue inhibitor of metalloproteases (Timp)

Product Code	CSB-MP894223DLU
Abbreviation	Timp
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VH14
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	CSC MPSHPQTHFA QADYVVQLRV LRKSDTIEPG RTTYKVHIKR TYKATSEARR MLRDGRLSTP QDDAMCGINL DLGKVYIVAG RMPTLNICSY YKEYTRMTIT ERHGFSGGYA KATNCTVTPC FGERCFKGRN YADTCKWSPF GKCETNYSAC MPHKVQTVNG VISRCRWRT QLYRKCMSNP
Source	Mammalian cell
Target Names	Timp
Protein Names	Recommended name: Tissue inhibitor of metalloproteases
Expression Region	28-210
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.