



Recombinant Drosophila melanogaster DNA-directed RNA polymerase II subunit RPB11 (Rpb11)

Product Code	CSB-MP894234DLU
Abbreviation	Rpb11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VJE4
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MNAPPTFESF LLYEGEKKII KELDTKVTNA AIFTINKEDH TLGNMIRNQL LKDPNVLFAG YKVPHPLEHK FVIRIQTAD YSPQEAFMNA ITDLLAELSL FEERFKDAIK EKKEGGD
Source	Mammalian cell
Target Names	Rpb11
Protein Names	Recommended name: DNA-directed RNA polymerase II subunit RPB11 Short name= RNA polymerase II subunit B11 Alternative name(s): DNA-directed RNA polymerase II 13.3 kDa polypeptide DNA-directed RNA polymerase II subunit J
Expression Region	1-117
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.