



Recombinant Drosophila melanogaster Protein Dr1 (NC2beta)

Product Code	CSB-EP894235DLU
Abbreviation	NC2beta
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VJQ5
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSNPQEELLP PSAEDDELTL PRASINKIIK ELVPTVRVAN ESRELILNCC SEFIHLISSE ANEVCNMRNK KTINAEHVLE ALERLGFHDY KQEAEAVLHD CKEVAAKRRR QSTRLENLGI PEEELLRQQQ ELFAKAREEQ AREEQQQWMS MQAAAMVQRP PLADGSVASK PSEDDDDDDD DDY
Source	E.coli
Target Names	NC2beta
Protein Names	Recommended name: Protein Dr1 Alternative name(s): Down-regulator of transcription 1 Negative cofactor 2-beta Short name= NC2-beta dNC2
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.