



Recombinant Drosophila melanogaster Cytosolic Fe-S cluster assembly factor NUBP2 homolog (CG4858)

Product Code	CSB-MP894259DLU
Abbreviation	CG4858
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VPD2
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MLDKVKNVIV VLSGKGGVGK STVSTQLSLA LRKNGFKVGL LDIDLCPSPV PYLLGLEGRD IFQCDDGWVP VYTDESQTLA VMSIGFLLKN REDPVIWRGP KKTMMIRQFL TDVRWDELDT LIIDTPPGTS DEHITVMECL KEVGCHGAIL VTTTPQEVALD DVRKEITFCK KTGINILGIV ENMSGFVCPH CTCTNIFSS NGGVSLATYA QVPHLGTLPD DPRVGILAGT TTSVLDELDP STTAEVLTHI VEKLKTMMLVS
Source	Mammalian cell
Target Names	CG4858
Protein Names	Recommended name: Cytosolic Fe-S cluster assembly factor NUBP2 homolog
Expression Region	1-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.