



Recombinant Drosophila melanogaster Alpha-tubulin N-acetyltransferase 1 (CG3967)

Product Code	CSB-BP894269DLU
Abbreviation	CG3967
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VSY4
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MVEFRFDIKP LFAQPIIKVT SNLLPNTFRG DRRQCLDATS KMTEIIDQLG QLSATSQGLS KPVTTAQRRLR MSDNQTIYLL ADNEAGHNGA VLGLLKVGTK NLYLFDEAGK TRMVEQTPSI LDFYVHESRQ RAGLGKRLFQ TMLNEEQWTA RKCSVDRPSE KLLSFLSKHY GLKRIIPQAN NFVLYEGFFN DGESGNGGGN GHANGTPNGL HITNSPNTHL FGATYLGEDS NQRRGSQQQT TPNARLQQIT QISPSGRYGA KRPTCSMAEI IHAGNSKGGN GNGSAEANSNG NDNHDIPEIA EQLQRQSLAD LEANSYEPEP EVEPEPEPEP EPEPEPEVIT PPSPPPKSHT PTPPSVRSPE VAESIVGDNR RAPAFQLSKQ HTGMKNRSFG VGMAVMPSSK MEFDQMERED FGVVKINRPP GHEVTSPGQD NTDAMSTVSS GGGGLTDQGY YDLKFYHNKL W
Source	Baculovirus
Target Names	CG3967
Protein Names	Recommended name: Alpha-tubulin N-acetyltransferase 1 Short name= Alpha-TAT 1 Short name= TAT 1 EC= 2.3.1.108 Alternative name(s): Acetyltransferase mec-17 homolog 1
Expression Region	1-461
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.