



Recombinant Chicken Ataxin-3 (ATXN3)

Product Code	CSB-MP894318CH
Abbreviation	ATXN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W689
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MESIFHERQE GSLCAQHCLN NLLQGEYFSP VELSSIAQQL DEEERMMAE GGVSSEYRT FLQQPSVNMD DSGFFSIQVI SNALKVWGLE LILFNSPEYQ RLGIDPINEK SFICNYKEHW FTVRKLGKQW FNLNSLLMGP ELISDTYLAL FLAQLQQEGY SIFVVKGDLP DCEADQLLQM IRVQQVQRPK LIGEETAQSR DQRLPRSDVD QAIEVSHPFDTGTGMLDEDEE NFQRALALSR QEIDMEDEEA DLRRAIQLSM QGSRQSEFSN SLPQNASQPP HTSQTDSLSS EDLRRRRQAY FEKQQQLQQ QDLTLNLHDK PTINSSTLEA DPGGDMSEED MLQAAMNMSL ESARNHLSTE EKK
Source	Mammalian cell
Target Names	ATXN3
Protein Names	Recommended name: Ataxin-3 EC= 3.4.19.12 Alternative name(s): Machado-Joseph disease protein 1 homolog
Expression Region	1-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	Machado-Joseph disease, also known as spinocerebellar ataxia-3, is an autosomal dominant neurologic disorder. This protein contains (CAG) _n repeats in the coding region, and the expansion of these repeats from the normal 13-36 to 68-79 is one cause of Machado-Joseph disease. There is a negative correlation between the age of onset and CAG repeat numbers. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.