



Recombinant Mouse RAC-gamma serine/threonine-protein kinase (Akt3)

Product Code	CSB-MP894361MO
Abbreviation	Akt3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WUA6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	SDVTIVKEG WVQKRGEYIK NWRPRYFLK TDGSFIGYKE KPQDVDLPYP LNNFSVAKCQ LMKTERPKPN TFIIRCLQWT TVIERTFHVD TPEEREEWTE AIQAVADRLQ RQEEERMNCS PTSQIDNIGE EEMDASTTHH KRKTMNDFDY LKLLGKGTFG KVILVREKAS GKYYAMKILK KEVIIAKDEV AHTLTESRVL KNTRHPFLTS LKYSFQTKDR LCFVMEYVNG GELFFHLSRE RVFSEDRTRF YGAEIVSALD YLHSGKIVYR DLKLENLMLD KDGHIKITDF GLCKEGITDA ATMKTFCGTP EYLAPEVLED NDYGRAVDWW GLGVVMYEMM CGRLPFYNQD HEKLFELILM EDIKFPRTLS SDAKSLLSGL LIKDPNKRLG GGPDDAKEIM RHSFFSGVNW QDVYDKKLV PFKPQVTSET DTRYFDEEFT AQTITITPPE KYDDDGMDGM DNERRPHFPQ FSYSASGRE
Source	Mammalian cell
Target Names	Akt3
Protein Names	Recommended name: RAC-gamma serine/threonine-protein kinase EC= 2.7.11.1 Alternative name(s): Protein kinase Akt-3 Protein kinase B gamma Short name= PKB gamma RAC-PK-gamma
Expression Region	2-479
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a member of the AKT, also called PKB, serine/threonine protein kinase family. AKT kinases are known to be regulators of cell signaling in response to insulin and growth factors. They are involved in a wide variety of biological processes including cell proliferation, differentiation, apoptosis, tumorigenesis, as well as glycogen synthesis and glucose uptake. This kinase has been shown to be stimulated by platelet-derived growth factor (PDGF), insulin, and insulin-like growth factor 1 (IGF1). Alternatively splice transcript



variants encoding distinct isoforms have been described.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.