



# Recombinant Rat Proto-oncogene tyrosine-protein kinase Src (Src)

<b>Product Code</b>	CSB-BP894363RA
<b>Abbreviation</b>	Src
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9WUD9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GSNKS PKD ASQRRRSLEP AENVHGAGGA FPASQTSPSKP ASADGHRGPN AAFVPPAAAE PKLFGGFNSS DTVTSPQRAG PLAGGVTTFFV ALYDYESRTE TDL SFKKGER LQIVNNTTEGD WWLAHSLSTG QTGYIPSNYV APSDSIQAE WYFGKITRRE SERLLLNAEN PRGTFLVRES ETTKGAYCLS VSDFDNAKGL NVKHYKIRKL DSGGFYITSR TQFNSLQQLV AYYSKHADGL CHRLTTVCPT SKPQTQGLAK DAWEIPRESL RLEVKLGQGC FGEVWMTWN GTTRVAIKTL KPGTMSPEAF LQEAQVMKKL RHEKLVQLYA VVSEEPYIV TEYMNKGSLL DFLKGETGKY LRLPQLVDMS AQIASGMAYV ERMNYVHRDL RAANILVGEN LVCKVADFGL ARLIEDNEYT ARQGAKFPIK WTAPEAALYG RFTIKSDVWS FGILLTELTT KGRVPYPGMV NREVLQDQVER GYRMPCPEEC PESLHDLMCQ CWRKEPEERP TFEYLQAFLE DYFTSTEPQY QPGENL
<b>Source</b>	Baculovirus
<b>Target Names</b>	Src
<b>Protein Names</b>	Recommended name: Proto-oncogene tyrosine-protein kinase Src EC= 2.7.10.2 Alternative name(s): Proto-oncogene c-Src pp60c-src Short name= p60-Src
<b>Expression Region</b>	2-536
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene is highly similar to the v-src gene of Rous sarcoma virus. This proto-oncogene may play a role in the regulation of embryonic development and cell growth. This protein is a tyrosine-protein kinase whose activity can be inhibited by phosphorylation by c-SRC kinase. Mutations in this gene could be involved in the malignant progression of colon cancer. Two transcript variants encoding the same protein have been found for this gene.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.