



Recombinant Zea mays Caffeoyl-CoA O-methyltransferase 2 (CCOAOMT2)

Product Code	CSB-YP894508ZAX
Abbreviation	CCOAOMT2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XGD5
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	MATTATEATK TTAPAQEQQA NGNGNGEQKT RHSEVGHKSL LKSDDLQYI LDTSVYPREP ESMKELREIT AKHPWNLMTT SADEGQFLNM LIKLIGAKKT MEIGVYTGYS LLATALALPE DGTILAMDIN RENYELGLPC INKAGVGHKI DFREGPALPV LDDLVDKEQ HGSFDFAFVD ADKDNYLNYH ERLKLVPRG GLIGYDNTLW NGSVLPDDA PMRKYIRFYR DFVLALNSAL AADDRVEICQ LPVGDGVTLG RRVK
Source	Yeast
Target Names	CCOAOMT2
Protein Names	Recommended name: Caffeoyl-CoA O-methyltransferase 2 EC= 2.1.1.104 Alternative name(s): Trans-caffeoyl-CoA 3-O-methyltransferase 2 Short name= CCoAMT-2 Short name= CCoAOMT-2
Expression Region	1-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.