



Recombinant Arabidopsis thaliana B3 domain-containing protein REM19 (REM19)

Product Code	CSB-YP894519DOA
Abbreviation	REM19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XIB5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MQMDSAQNQF NKRARLFEDP ELKDAKVIYP SNPESTEPVN KGYGGSTAIQ SFFKESKAEE TPKVLKKRGR KKKNPNPEEV NSSTPGGDDS ENRSKFYESA SARKRTVTAE ERERAVNAAK TFEPTNPYFR VVLRPSYLYR GCIMYLPSGF AEKYLSGISG FIKLQLGEKQ WPVRCLYKAG RAKFSQGWYE FTLENNIGEG DVCVFELLRT RDFVLEVTAF RVNEYV
Source	Yeast
Target Names	REM19
Protein Names	Recommended name: B3 domain-containing protein REM19 Alternative name(s): Protein RELATED TO VERNALIZATION 1 Protein REPRODUCTIVE MERISTEM 19
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.