



Recombinant Probable cytosolic iron-sulfur protein assembly protein CIAO1 homolog (Y18D10A.9)

Product Code	CSB-YP894564CXY
Abbreviation	Y18D10A.9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XW12
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MLRQIGEFYH QGEKDDTSRV WMTCWHHGGR ILASCGDDKA VRVWSLVGEP DSKMRLECRT TLDDSHTRAV RSVAFSNDGK CLVSASF DAS VVYQQEDGE FAEVNKLEGH ESEVKCAVFS KSDEF LATCS RDKSVWFQQ DEDEDFSVSS ILQPHTQDVK QVAWHPTEDL LVSCSYDSSI RFYRFDGEDW VTQQKIDGCH VGTVWSIAFD TEGHRLVTVG EDHCIQLFVR ENIGSKSADQ DTWKS VARYD VENTRWPLYS VAWNSTNDVI ATGGGDCKIR LFKISSTPES PVIEHLGVVG RHELDVNHVA WNP NPKFSNL LTSASDDGTI RLWELEI
Source	Yeast
Target Names	Y18D10A.9
Protein Names	Recommended name: Probable cytosolic iron-sulfur protein assembly protein CIAO1 homolog
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.