



Recombinant Probable actin-related protein 2/3 complex subunit 3 (arx-5)

Product Code	CSB-MP894566CXY
Abbreviation	arx-5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XWV3
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MPAYHSKFDT ELKVLPLGNT NMGKLPVRTN FKGPAQTNQ DDIIDEALTY FKPNIFFREF EIKGPADRTM IYLIFYITEC LRKLQKSPNK IAGQKDLHAL ALSHLLPIPG ENGFLNSMY KAPQSKPDED EMRAYLQQIR QEIGARLCDL AFPDPQDRPS KWWLCFARRR FMDKGLVGQG VNL
Source	Mammalian cell
Target Names	arx-5
Protein Names	Recommended name: Probable actin-related protein 2/3 complex subunit 3 Alternative name(s): Arp2/3 complex 21 kDa subunit Short name= p21-ARC
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.