



Recombinant *Drosophila melanogaster* Protein CWC15 homolog (c12.1)

Product Code	CSB-YP894635DLU
Abbreviation	c12.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V3B6
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MTTAARPTFD PARGGSGRGE KDLSALSKQY SSRDLPGHTK LKYRETGQGT SDENRNRDFR KELEEREREA RSGTGATSSS SGKALPSIVR KAIEANNAGG GSSAAKRSKP DAGQQQAQQA AQQQAANMDA DEPLDNDSSD SDSDSDDDDA ALLAELQKIK QERLQETARR ESEKKQEDER IRMENILSGN PLMNYEPGTA ASAAGRASGL GGDLKIKRRW DDDVVFKNCA RSAPDKKTHF VNDALRSDFH KKFMDKYIK
Source	Yeast
Target Names	c12.1
Protein Names	Recommended name: Protein CWC15 homolog
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.