



# Recombinant *Drosophila melanogaster* 26S proteasome non-ATPase regulatory subunit 6 (Rpn7)

<b>Product Code</b>	CSB-BP894637DLU
<b>Abbreviation</b>	Rpn7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9V3G7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPAENLEEQG LEKNPNLELA QTKFLLTLAE YKQDAALKAK LLEAIRTENM APWYEHICSE LGWTVDKDLL ARMKENNRVE VEQLDAAIED AEKNLGEMEV REANLKKSEY LCRIGDKAAA ETAFRKTYEK TVSLGHRLDI VFHLIRLGLF YLDHDLITRN IDKAKYLIEE GGDWDRRNL KVYQGVSVA VRDFKAAATF FLDTVSTFTS YELMDYPTFV RYTVYVAMIA LPRNELRDKV IKGSEIQEVL HGLPDVKQFL FSLYNQCQYEN FYVHLAGVEK QLRLDYLIHP HYRYVREMR ILGYTQLLES YRSLTLQYMA ESFGVTVEYI DQELARFIAA GRLHAKVDRV GGIVETNRPD NKNWQYQATI KQGDLLLNI QKLSRVINI
<b>Source</b>	Baculovirus
<b>Target Names</b>	Rpn7
<b>Protein Names</b>	Recommended name: 26S proteasome non-ATPase regulatory subunit 6 Alternative name(s): 26S proteasome regulatory subunit S10 Protein Rpn7
<b>Expression Region</b>	1-389
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.