



Recombinant Drosophila melanogaster G patch domain and ankyrin repeat-containing protein 1 homolog (CG8152)

Product Code	CSB-MP894653DLU
Abbreviation	CG8152
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V7A7
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MNSQNELHPN WLALTTLPLQ LKRFVRAGDP DFIAESPKAV KHQIDGLNGA EAREFYKEVL DAP TTCQPNL PSTHKT AQTR REPERILVPF DKSKFFRLAT SNKVEELSQM KASEEDLNSR DSGW TALMM ACEGATEAV SWLVQRGVQV ETS DKSGNTA LKLAQRKGHL DVVHLLLESLP ILEETSEEDE SVDGN NPFYC EICKRDYKET PWPIHQ TSTV HQFN LKALPA HKLHKFNISA KNRGLQLMVK QGWDQEHLG PSQSGRLYPV KTVLRKQRTG LGIEQQSARV SHFGAFDLNA VRRRDPIYQP RRTRS DMQRE KVREWKRERY LRRVLS
Source	Mammalian cell
Target Names	CG8152
Protein Names	Recommended name: G patch domain and ankyrin repeat-containing protein 1 homolog
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.