



Recombinant Drosophila melanogaster Ragulator complex protein LAMTOR2 homolog (CG5189)

Product Code	CSB-EP894655DLU
Abbreviation	CG5189
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V8I2
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MLKPKALTQV LSQANTGGVE NTLLLSQEGA LLAYSGYGDK DARITAAIAS NIWAAYEKHG RNAFREGRLT FVLIDCENGH VAITQVASVL LCLYAKQTVG LGLLKQKAMS LASYLERPLK QISAS
Source	E.coli
Target Names	CG5189
Protein Names	Recommended name: Ragulator complex protein LAMTOR2 homolog
Expression Region	1-125
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.