



# Recombinant Drosophila melanogaster Interleukin enhancer-binding factor 2 homolog (CG5641)

<b>Product Code</b>	CSB-MP894688DLU
<b>Abbreviation</b>	CG5641
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9VG73
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Drosophila melanogaster (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVRGALRGGR PMRGGIRPPF KKTFVPRHPF DLTAEVFFP KVPSAGAVDD SALTAALLKR NQDLSPTPSE QTAIGNLVTK VQAVLDNLVV APGDLTTCQL EEVRQVGSFK KGTILTGNV ADVVVILKTL PTKEAVDALA KKVEADLKAS MKTEVLTKGD QHTVQIHERG FDIANVHAKV RILIATLPQN LRKLEPEIHL DHKLMQSHLA AIRHTRWFEE NAHHSSIKVL IRILKDLTRR FDAFSPLSAW MLDLIAHLAI MNNPSRQALP INLAFRRVFQ LLSAGLFLPG SAGITDPTPEP GHIRVHTAMT LEQQDVCCYT SQTLLRVLAH GGYKHILGLE GNTSVVREMS VWNGVCISPL TAVYEKPTDK KEGDLEEDID MIENENEEEG SDDGAE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CG5641
<b>Protein Names</b>	Recommended name: Interleukin enhancer-binding factor 2 homolog
<b>Expression Region</b>	1-396
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.