



Recombinant *Drosophila melanogaster* Glutathione S-transferase D4 (GstD4)

Product Code	CSB-YP894689DLU
Abbreviation	GstD4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VG96
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDFYSPRSS GSRTIIMVAK ALGLELNKKQ LRITEGEHLK PEFLKLN PQH TIPTLVDNGF AIWESRAIAV YLVEKYGKDD SLFPNDPQKR ALINQRLYFD MGTLDHDSFMK YYYPFIRTGQ LGNAENYKKV EAAFEFLDIF LEGQDYVAGS QLTVADIAIL SSVSTFEVVE FDISKYPNVA RWYANAKKIT PGWDENWKGL LQMKTMYEAQ KASLK
Source	Yeast
Target Names	GstD4
Protein Names	Recommended name: Glutathione S-transferase D4 Short name= DmGST23 EC= 2.5.1.18
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.