



# Recombinant *Drosophila melanogaster* Elongation factor Ts, mitochondrial (mEFTs)

<b>Product Code</b>	CSB-EP894701DLU
<b>Abbreviation</b>	CG6412
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9VJC7
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AG SGGSGGLEKS ALAALRKKTG YTFANCKKAL EKHNNNDVGLA EKWLHEQAQT LGWSKATKVA DRATAHGLIG VLIRGNRGAM VELNCETDFV ARNDTFKRFV DHVACMCLQY TDLTDFDGD L WKLGFADADAL RNLRTTEGRT LGDHLALLIG AIGENATIRR ALCFKANNDL KLVGYAHPAP TNVGTTEGIT QVGKYGAIVA YRSTHPLLDF EFHKSICQQI VGMKPTKIGE YDKDKPAENK DDETCLIHQE YLLDADKTVG EALQEHNCEI VDYHRFECGE HTERSLEAII RSQHQSSN
<b>Source</b>	E.coli
<b>Target Names</b>	mEFTs
<b>Protein Names</b>	Recommended name: Elongation factor Ts, mitochondrial Short name= EF-Ts Short name= EF-TsMt
<b>Expression Region</b>	19-318
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.