



Recombinant Drosophila melanogaster Cytosolic Fe-S cluster assembly factor NUBP1 homolog (CG17904)

Product Code	CSB-MP894702DLU
Abbreviation	CG17904
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VJI9
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MQAPPPEHCP GVESEEAGKG SACSGCPNQG LCSDPNKKLE DPGKALVVES MKDVKHKLLI LSGKGGVGKS TVTSLTRYL ARSNPDSNFG VLDIDICGPS QPRLMGALGE SVHQSGYGWS PVGIEDNVCL MSIGFLLGSV DDAIWRGPK KNGMIRQFLS EVDWGNLDLL LLDTPPGTSD EHLSVVSYLK DDANPESLRA VMVTTPQEV S LLDVRKEINF CKKQNIPIVG VIENMSSFRC GHCGNSSEIF PAKTGGAAAM CAEMGIPLLG SLPLDQQISK ACDSGEDLTE FKNVTTEALE GICSKIMASF S
Source	Mammalian cell
Target Names	CG17904
Protein Names	Recommended name: Cytosolic Fe-S cluster assembly factor NUBP1 homolog
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.