



# Recombinant *Drosophila melanogaster* Selenide, water dikinase 2 (Sps2)

<b>Product Code</b>	CSB-YP894706DLU
<b>Abbreviation</b>	Sps2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9VKY8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MFQPEKHGLE PDFQLTKFTT HTGUSCKIPQ KVLEKYLRGT EIENKNNDGY LIGSGMDCAV IPLKRHKDYL LIQTVDFYYP MVNDPELLGR IALANVLSDV YAVGVTQFDT VEMIVSTSTS FSEKERDVVI GLVMKGFQNS LKANGYRNTNTP LIIRQLKINP WCIIGGIATS VCRSEEIILP SNAQPGDVLV LTKPLGGQMA MDAHLWQLNQ TEKYKLLSE CSDADIKETF EIAVKSMTYL NKNAALLMHK YQAHCATDIT GFGLLGHANN LAQFQKEKVL FQINKLPIIK NVLKFSTLVG QSTKFRSGRS VETSGGLLIC LPADAADKFC RDFEETNGE QKSFQIGHVT AANESDAVLC EDVEFIEVSL
<b>Source</b>	Yeast
<b>Target Names</b>	Sps2
<b>Protein Names</b>	Recommended name: Selenide, water dikinase 2 EC= 2.7.9.3 Alternative name(s): Selenium donor protein 2 Short name= Dsps2 Selenophosphate synthase 2
<b>Expression Region</b>	1-370
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.