



Recombinant *Drosophila melanogaster* Probable queuine tRNA-ribosyltransferase (Tgt)

Product Code	CSB-MP894721DLU
Abbreviation	Tgt
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VPY8
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MGPSHIPPLT YKVVAECSVS KARAGLMTLR HSEVNTPVFM PVGTQGTLKG IVPDQLIELN CQILLGNTYH LGLRPGIETL KKAGGLHKFM GWPRAILTDS GGFQMVSLQ LAEIDEHGVN FRSPFDNSQC MLTPEHSIEI QNAIGGDIMM QLDDVVKTTT TGPRVEEAME RTIRWVDRCI EAHARDDDQS LFPIVQGGGLD VPLRQRCVSA LMERQVRGFA VGGLSGGESK HDFWRMVDVC TGYPKDKPR YLMGVGFAAD LVVCVALGID MFDCVFPTRT ARFGCALVDS GQLNLKQPKY KLDMEPIDKD CDCSTCRRYT RSYLHHIATN ESVSSLLSI HNVAYQLRLM RSMREAIQRD EFPQFVADFM ARHFKAEPVP AWIREALSAV NIQLPADPER IDEQDQPKT EKRRETEDVA EEQVASS
Source	Mammalian cell
Target Names	Tgt
Protein Names	Recommended name: Probable queuine tRNA-ribosyltransferase EC= 2.4.2.29 Alternative name(s): Guanine insertion enzyme tRNA-guanine transglycosylase
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.