



Recombinant *Drosophila melanogaster* Nuclear distribution protein nudE homolog (nudE)

Product Code	CSB-EP894734DLU-B
Abbreviation	nudE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VT70
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MESPPMFNSV EDECRYWKER SKQYHKEWTD VKQEYDEFVE QSREMEIEMD ATLDQKQSII KDLTAKLTMF ERENESLKLK LESHGIDMSN MEKQLETVKK DRDTMKVYLR QLEQKNDDLE RAHRILNESI ENFEKMLDQA YEKNALLELE VDEKGLLQEK LQRLMDETRD LKQELNVKSR FTPVVNGTSV PTANDTNTVN SSMNSSASLP NGIVANGELV KHDNAVATRA TSVSVNALNG SLVNRNEYNQ QHSLKSDGDE TEDVPREWSA NAAATPLPA IHTSGAGRTA LLLAGCRVWV ILQFSVIHKP QIRNPIT
Source	<i>E.coli</i>
Target Names	nudE
Protein Names	Recommended name: Nuclear distribution protein nudE homolog
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.