



Recombinant *Drosophila melanogaster* Accessory gland protein Acp76A (Acp76A)

Product Code	CSB-EP894744DLU-B
Abbreviation	Acp76A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VWW1
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	QNVSFQLI REIDRYTPEN FVLSVLNIEM ILFEIHAACA VESNNDLERS LIINFGYSEA RQEVLDWGLR YKKASSAKFQ MANKVAVSQK LPLSQKLRLLV NEVLM TSAKK YDVTKDVRPS KLMDEWLSSH LDGVLANFVQ EKKLNAGENI VAISGMTVTP LWASHFQSEI NRYFVNPNPGT GYASKDPTCV PMMHSLSSFE TMSTDEAKGI YIPFSSANLG MLILLPRKGV TCKDILDNLN NQINVEYNDH KDVHLLLPIF KEKFDYNIK FFNGINIEDT FKDSAFKSKA KIKINFRVN HGIRFQPILR LEVDDIDTG KTETFEVNRPFVFKDKIN VYAVGRIENL DGLTDKVNCS KKYADL
Source	E.coli
Target Names	Acp76A
Protein Names	Recommended name: Accessory gland protein Acp76A
Expression Region	23-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.