



Recombinant *Cricetulus griseus* NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 1 (NDUFA1), partial

Product Code	CSB-BP894818DXU
Abbreviation	NDUFA1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WU08
Product Type	Recombinant Protein
Immunogen Species	<i>Cricetulus griseus</i> (Chinese hamster) (<i>Cricetulus barabensis griseus</i>)
Purity	>85% (SDS-PAGE)
Source	Baculovirus
Target Names	NDUFA1
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 1 Alternative name(s): Complex I-MWFE Short name= CI-MWFE NADH-ubiquinone oxidoreductase MWFE subunit
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.