



Recombinant Mouse 6-phosphofructokinase type C (Pfkp)

Product Code	CSB-MP894822MO
Abbreviation	Pfkp
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WUA3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSDLDS SSSSS AYPKYLEHLS GDGKAIGVLT SGGDAQGMNA AVRAVVRMGI YTGAKVYFIY EGYQGLVDGG SNIVEAKWDC VSSILQVGGT IIGSARCKAF RSREGRLKAA CNLARLGITN LCVIGGDGSL TGANLFRKEW SGLLEELARN GDIDNDTVQK YSYLNVVGMV GSIDNDFCGT DMTIGTDSAL HRIIEVVDI MTTAQSHQRT FVLEVMGRHC GYLALVSALT CGADWVFLPE SPPEEDWEEN MCLKLSENRA RKKRLNIIIV SEGAIMQNK PITSEKIKEL VVKNLGFDR VTILGHVQRG GTPSAFDRIL ASRMGVEAVI ALLEATPETP ACVVSLRGNQ AVRLPLMECV QMTQDVQKAM DERRFKEAVK LRGRRFEGNL NTYKRLAIKL PDEKIVKSNC NVAVINVGAP AAGMNAAVRS AVRVG IADGH KMFAIYDGE GFANGQIKEI GWADVGGWTG QGGSILGTKR TLP GKYLEKI AEQMHSHSIN ALLIIGGFEA YLGLLELAAA REKHEAFCVP MVMVPATVSN NVP GSDFSIG ADTALNTITD TCDRIKQSAS GTKRRVFIIE TMGGYCYLA NMGALAAGAD AAYIFEEPFD IGD LQSNVH L TEKMKTSIQ RGLVLRNESC SVNYTTDFIY QLYSEEGKGV FDCRKNVLGH MQQGGAPSPF DRNFGTKISA KAMEWISAKL KGSQGTGKKF VSDDSICVLG ICKRDLLFQP VAELKKVTDF EHRIPKEQWW LKL RPIMKIL AKYEASYDMS DSGKLESLQH HEEL
Source	Mammalian cell
Target Names	Pfkp
Protein Names	Recommended name: 6-phosphofructokinase type C EC= 2.7.1.11 Alternative name(s): Phosphofructo-1-kinase isozyme C Short name= PFK-C Phosphofructokinase 1 Phosphohexokinase
Expression Region	1-784
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.