



Recombinant Arabidopsis thaliana 3-phosphoinositide-dependent protein kinase 1 (PDPK1)

| | |
|--------------------------|---|
| Product Code | CSB-MP894967DOA |
| Abbreviation | PDPK1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9XF67 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MLAMEKEFDS KLVLQGNSSN GANVSRSKSF SFKAPQENFT SHDFEFGKIY GVGSYSKVVR AKKKETGTVY ALKIMDKKFI TKENKTAYVK LERIVLDQLE HPGIKLYFT FQDTSSLYMA LESCEGGELF DQITRKGRLS EDEARFYTAE VVDALYIHS MGLIHRDIKP ENLLTSDGH IKIADFGSVK PMQDSQITVL PNAASDDKAC TFGTAAYVP PEVLNSSPAT FGNDLWALGC TLYQMLSGTS PFKDASEWLI FQRIIARDIK FPNHFSEAR DLIDRLDTE PSRRPGAGSE GYVALKRHPF FNGVDWKNLR SQTPPKLAPD PASQTASPER DDTHGSPWNL THIGDSLATQ NEGHSAPPTS SESSGSITRL ASIDSFDSRW QQFLEPGESV LMISAVKKLQ KITSKKVQLI LTNKPKLIYV DPSKLVVKGK IIVSDNSNDL NVVVTSPSHF KICTPKKVLS FEDAKQRASV WKKAJETLQN R |
| Source | Mammalian cell |
| Target Names | PDPK1 |
| Protein Names | Recommended name: 3-phosphoinositide-dependent protein kinase 1 Short name= AtPDK1 EC= 2.7.11.1 |
| Expression Region | 1-491 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.