



# Recombinant Horse Chromogranin-A (CHGA)

<b>Product Code</b>	CSB-YP895006HO
<b>Abbreviation</b>	CHGA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9XS63
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Equus caballus (Horse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	LP VNSPMDTGDT EVMKCIVEVI SDTLSKPSPV PVSQECFETL RGDERILSIL RHQNLKELQ DLALQGAKER APQKHSRLE DELAEVLEKQ NHQAEKKEVT EEALSEDAAE ARGDSKEVEE NGEDADGARP QAALEPEQES RVEDAQAPGE EKEAINTHSP TRLP SQKHPD PQAEGDSDSP SQGLVDREKG LGAERGQQA REEEEDAEAGE KADAE EEGPT AAFNPHPSLS YKIRKGESWS EALVVDGARK TGAEAAQPPE GQGEREHSRQ EEEEEETAG ASRGLFRGGK SRELEQEKEQ ERLSKEWEDA KRWSKMDQLA KELTAEKRLE GEDEEEDDPD RSMKLSFRAR AYGFRG PGLQ LRRGWRPSSR EDSIEAGLPP PVRGYPEEK EEEGSANRRP EDQELESLSA IEAELEKVAH QLQALRRG
<b>Source</b>	Yeast
<b>Target Names</b>	CHGA
<b>Protein Names</b>	Recommended name: Chromogranin-A Short name= CgA Cleaved into the following 2 chains: 1. Pancreastatin 2. WE-14
<b>Expression Region</b>	19-448
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. It is found in secretory vesicles of neurons and endocrine cells. This gene product is a precursor to three biologically active peptides; vasostatin, pancreastatin, and parastatin. These peptides act as autocrine or paracrine negative modulators of the neuroendocrine system. Other peptides, including chromostatin, beta-granin, WE-14 and GE-25, are also derived from the full-length protein. However, biological activities for these molecules have not been shown.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.