



Recombinant *Saimiri boliviensis boliviensis* Peptidyl-prolyl cis-trans isomerase FKBP5 (FKBP5)

Product Code	CSB-MP895009SWC
Abbreviation	FKBP5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XSH5
Product Type	Recombinant Protein
Immunogen Species	<i>Saimiri boliviensis boliviensis</i> (Bolivian squirrel monkey)
Purity	>85% (SDS-PAGE)
Sequence	MTTDEGAKNS RGNPAATVAE QGEDVTSKGD RGVLKIVKRV GHGEETPMIG DRVYVHYNGK LANGKKFDSS HDRNEPFVFS IGKGQVIKAW DIGVATMKKG EICHLCKPE YAYGATGSLP KIPSNATLFF EVELLDFKGE DLLEDGGIIR RTKRRGEGYS NPNEGARVQI HLEGRCGGRV FDCRDVAFTV GEGEDHDIPI GIDKALEKMQ REEQCILHLG PRYGFGEAGK PKFGIEPNAE LIYEVTLKSF EKAKESWEMD TKEKLEQAAI VKEKGTVYFK GGKYVQAVIQ YGKIVSWLEM EYGLSEKESK ASESFLLAAF LNLAMCYLKL REYTKAVECC DKALGLDSAN EKGLYRRGEA QLLMNEFESA KGDFEKVLEV NPQNKAARLQ IFMCQKKAKE HNERDRRTYA NMFKKFAEQD AKEEANKAMS KKTSEGVTNE KLTASHAVEE EKPEGHV
Source	Mammalian cell
Target Names	FKBP5
Protein Names	Recommended name: Peptidyl-prolyl cis-trans isomerase FKBP5 Short name= PPlase FKBP5 EC= 5.2.1.8 Alternative name(s): 51 kDa FK506-binding protein Short name= 51 kDa FKBP Short name= FKBP-51 FK506-binding protein 5
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.