



# Recombinant Monodelphis domestica L-lactate dehydrogenase A chain (LDHA)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP895018MUE   |
| <b>Abbreviation</b>      | LDHA  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9XT87  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Monodelphis domestica (Gray short-tailed opossum)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | GTVKDQLIL NVLKEEQTPH NKITVVGVA VGMACAISIL MKDLADELAL<br>VDVIEDKLLG EMMDLQHGS FLKTPKIVSS KDYAVTANSK LVVITAGARQ<br>QEGESRLNLV QRNVNIFKFI IPNIVKYSPN CKLLVVSNPV DILTYVAWKL<br>SGFPKNRVIG SGCNLD SARF RYLMGEKLG I HSSSCHGWIL GEHGDSSVPV<br>WSGVNVAGVS LKSLHPALGT DSDSEQWKDV HKQVVESAYE VIKLKGYSW<br>AIGLSVADLA ESIMKNLRRV HPISTMIKGL YGINEDVFLS VPCILGQNGI<br>SDVVKVTLTT EEESRLKQSA DTLWGIQKEL QF               |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | LDHA  |
| <b>Protein Names</b>     | Recommended name: L-lactate dehydrogenase A chain Short name= LDH-A<br>EC= 1.1.1.27 Alternative name(s): LDH muscle subunit Short name= LDH-M   |
| <b>Expression Region</b> | 2-332   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |