



Recombinant BTB and MATH domain-containing protein 45 (bath-45)

Product Code	CSB-BP895025CXY
Abbreviation	bath-45
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XV51
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MAAAHKVFEL SHVFKDVSKI GDESVCSN EDYFGVPWKI VLRHKDEKMG VFLSCEKPSN NDAHCVNVV YSVKLMSSSG VCIKSDSKL FGEITEHGWV NFVNWNTMVD AFLIDDSIII EVLVKTIFMS GLPKKSPKNF DESNKEFSDG IVVVKHQNFH ILKKFLAFHC EYLEKLFFGD FKEAGKAEVT LKAISSTDFH YLLAVLYEDY AIDDDNVGGI LRLASYLQVP IVIRKCEQFL IEGSGKPMEI RLDIAEKECL QNLKKHCLSK VGTADQIEAA RAILPSDIPK RALKKFDESN KEFSDVIFVV EHNKFYVLKQ ILASHSSHFK KIFVENVDKR EFTLADIDSN DFQNLVEVMY GVSFIDDLTV EGVHLAHMY DRPFPMQKCV EFLIDLSEKS PKEKLIKA YQLKNLEKAA LARINLPESI RAALSCDMAQ MEVEVLKALM QKALFHLSDT KVHTETFF
Source	Baculovirus
Target Names	bath-45
Protein Names	Recommended name: BTB and MATH domain-containing protein 45
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.